

Anti-Cardiac Troponin I Monoclonal Antibody, Capture

AGMB06041

Description

This Anti-Cardiac Troponin I Monoclonal Antibody, Capture is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06041
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	TNNI3, TNNC1, Troponin I, cardiac muscle, Cardiac troponin I
Applications:	ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	-
Clonality:	Monoclonal Antibody
SwissProt ID:	P19429
Immunogen:	-
Gene ID:	7137
Gene Name:	TNNI3

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Mouse
Isotype:	-
Purification:	Affinity Purified
Conjugated:	-
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in PBS containing 0.09% sodium azide
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.