

# Anti-TARC/CCL17 Monoclonal Antibody,Capture (PBS Only)

AGMB06043

## Description

---

This Anti-TARC/CCL17 Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB06043
<b>Contents:</b>	100µg Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	TARC, CCL17, SCYA17, C-C motif chemokine 17, C-C motif chemokine 17, Small-inducible cytokine A17, Thymus and activation-regulated chemokine
<b>Applications:</b>	<b>ELISA</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	-
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q92583
<b>Immunogen:</b>	-
<b>Gene ID:</b>	6361

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** CCL17  
**Host Species:** Rabbit  
**Isotype:** -  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** SStore at 4°C.Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 0.01M PBS, pH 7.2

**Antibody Dilution Ratio:** Optimize antibody dilution ratio according to assay requirements and sample type.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.