

# Anti-BOB1 [4G9-4A8-3E8] Monoclonal Antibody

## AGMB06058

### Description

---

This Anti-BOB1 [4G9-4A8-3E8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB06058
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	OBF1, OCAB, OBF-1
<b>Applications:</b>	<b>IHC-P</b> <b>ELISA</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse,
<b>Clone:</b>	4G9-4A8-3E8
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q16633
<b>Immunogen:</b>	Synthesized peptide derived from human BOB.1
<b>Gene ID:</b>	5450
<b>Gene Name:</b>	POU2AF1

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse

**Isotype:** IgG2a

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:200
	ELISA	1:500-1:5000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.