

Anti-Cytokeratin 5/6 [5E7-10D10-2C6] Monoclonal Antibody

AGMB06060

Description

This Anti-Cytokeratin 5/6 [5E7-10D10-2C6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06060
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	58 kDa cytokeratin, Cytokeratin-5, CK-5, Keratin-5, K5, Type-II keratin Kb5, Cytokeratin-6A, CK-6A, Cytokeratin-6D, CK-6D, Keratin-6A, K6A, Type-II keratin Kb6, Hom s 5
Applications:	IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	5E7-10D10-2C6
Clonality:	Monoclonal Antibody
SwissProt ID:	P13647/P02538
Immunogen:	Synthesized peptide derived from human Cytokeratin 5/6
Gene ID:	3852/3853

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: KRT5/KRT6A

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.