

# Anti-CD105 [6A6-5B7-8G9] Monoclonal Antibody

AGMB06086

## Description

---

This Anti-CD105 [6A6-5B7-8G9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB06086
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	ENG, END, ORW, HHT1, ORW1, CD105, FLJ41744, Endoglin, CD antigen CD105, CD105/Endoglin
<b>Applications:</b>	<b>IHC-P</b> <b>ELISA</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	6A6-5B7-8G9
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P17813
<b>Immunogen:</b>	Synthesized peptide derived from human CD105(Endoglin)
<b>Gene ID:</b>	2022

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** ENG  
**Host Species:** Mouse  
**Isotype:** IgG2a  
**Purification:** Affinity Chromatography  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

	Application	Antibody Dilution Ratio
<b>Antibody Dilution Ratio:</b>	IHC	1:50-1:200
	ELISA	1:500-1:5000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.