

Anti-Cytokeratin 15 [2G8-10H10-8C9] Monoclonal Antibody

AGMB06087

Description

This Anti-Cytokeratin 15 [2G8-10H10-8C9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06087
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CK-15, CK15, cytokeratin15, K15, K1C15, K1CO, keratin 15, Keratin-15, keratin-15, beta, KRT15, KRTB, type I cytoskeletal 15
Applications:	IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	2G8-10H10-8C9
Clonality:	Monoclonal Antibody
SwissProt ID:	P19012
Immunogen:	Synthesized peptide derived from human Cytokeratin 15
Gene ID:	3866

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: KRT15
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.