

Anti-HER2 [3B4-9E2-7E4] Monoclonal Antibody

AGMB06102

Description

This Anti-HER2 [3B4-9E2-7E4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06102
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	NEU, NGL, ERBB2, TKR1, CD340, HER-2, VSCN2, MLN 19, MLN-19, c-ERB2, c-ERB-2, HER-2/neu, p185(erbB2)
Applications:	IHC-P
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	3B4-9E2-7E4
Clonality:	Monoclonal Antibody
SwissProt ID:	P04626
Immunogen:	Synthesized peptide derived from human Her-2(c-erbB-2)
Gene ID:	2064

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ERBB2

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.