

Anti-GCDFP 15 [4C6-8D2-1H1] Monoclonal Antibody

AGMB06115

Description

This Anti-GCDFP 15 [4C6-8D2-1H1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06115
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	BRST2, GCDGP15, gp17, GPIP4, PIP, Prolactin induced protein, SABP, Secretory actin binding protein
Applications:	IHC-P
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	4C6-8D2-1H1
Clonality:	Monoclonal Antibody
SwissProt ID:	P12273
Immunogen:	Synthesized peptide derived from human Gross cystic disease fluid protein-15 (GCDFP-15)
Gene ID:	5304

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: PIP

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.