

Anti-Mucin 16 [7A8-6E10-1H8] Monoclonal Antibody

AGMB06137

Description

This Anti-Mucin 16 [7A8-6E10-1H8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06137
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MUC16, CA125, Mucin-16, MUC-16, Ovarian cancer-related tumor marker CA125, CA-125, Ovarian carcinoma antigen CA125
Applications:	IHC-P
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	7A8-6E10-1H8
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8WXI7
Immunogen:	Synthesized peptide derived from human CA-125
Gene ID:	94025

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: MUC16

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.