

Anti-Progesterone Receptor [5C10-6H2-6B1] Monoclonal Antibody

AGMB06149

Description

This Anti-Progesterone Receptor [5C10-6H2-6B1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06149
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	PGR, NR3C3, Progesterone receptor, PR, Nuclear receptor subfamily 3 group C member 3
Applications:	IHC-P ICC/IF ELISA
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	5C10-6H2-6B1
Clonality:	Monoclonal Antibody
SwissProt ID:	P06401
Immunogen:	Synthesized peptide derived from human Progesterone Receptor(PR)
Gene ID:	5241

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: PGR
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:200
	IF	1:50-1:200
	ELISA	1:500-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.