

Anti-Bid [ZN139] Monoclonal Antibody

AGMB06200

Description

This Anti-Bid [ZN139] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06200
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FP497, MGC15319, MGC42355, BID, Apoptic death agonist, Apoptotic death agonist BID, BH3 interacting domain death agonist,
Applications:	IHC-P
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	ZN139
Clonality:	Monoclonal Antibody
SwissProt ID:	P55957
Immunogen:	Full length human recombinant protein of human BID(NP_001187) produced in HEK293T cell.
Gene ID:	637

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: BID

Host Species: Mouse

Isotype: IgG2b

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.