

Anti-CD8 alpha [ZN241] Monoclonal Antibody

AGMB06207

Description

This Anti-CD8 alpha [ZN241] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06207
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD8A, MAL, T-cell surface glycoprotein CD8 alpha chain, T-lymphocyte differentiation antigen T8/Leu-2, CD antigen CD8a
Applications:	IHC-P ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	ZN241
Clonality:	Monoclonal Antibody
SwissProt ID:	P01732
Immunogen:	Human recombinant protein fragment of human CD8A (NP_001759) produced in E.coli.
Gene ID:	925

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CD8A
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:200-1:400
	ELISA	1:500-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.