

Anti-Nkx3.1 [ZN195] Monoclonal Antibody

AGMB06253

Description

This Anti-Nkx3.1 [ZN195] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06253
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	NKX3.1, NKX3-1, NKX3.1, NKX3A, Homeobox protein Nkx-3.1, Homeobox protein NK-3 homolog A
Applications:	IHC-P
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	ZN195
Clonality:	Monoclonal Antibody
SwissProt ID:	Q99801
Immunogen:	Full length human recombinant protein of human NKX3-1 (NP_006158) produced in E.coli.
Gene ID:	4824

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: NKX3-1
Host Species: Mouse
Isotype: IgG2a
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.