

Anti-SERPINB4 [ZN16] Monoclonal Antibody

AGMB06262

Description

This Anti-SERPINB4 [ZN16] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06262
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	LEUPIN, PI11, SCCA-2, SCCA1, SCCA2, Serpin B4
Applications:	IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	ZN16
Clonality:	Monoclonal Antibody
SwissProt ID:	P48594
Immunogen:	Full length human recombinant protein of human SERPINB4 (NP_002965) produced in HEK293T cell.
Gene ID:	6318

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: SERPINB4

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.