

# Anti-VILIP [ZN115] Monoclonal Antibody

AGMB06272

## Description

---

This Anti-VILIP [ZN115] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB06272
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	Visinin-like 1, HLP3, VILIP, HPCAL3, HUVISL1, VILIP-1
<b>Applications:</b>	<b>IHC-P</b>
<b>Research Area:</b>	Signal Transduction
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	ZN115
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P62760
<b>Immunogen:</b>	Human recombinant protein fragment of human VSNL1 (NP_003376) produced in E.coli.
<b>Gene ID:</b>	7447

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** VSNL1  
**Host Species:** Mouse  
**Isotype:** IgG1  
**Purification:** Affinity Chromatography  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

<b>Antibody</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
<b>Dilution Ratio:</b>	IHC	1:50-1:100

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.