

# Anti-FOXP3 [ZN248] Monoclonal Antibody

## AGMB06293

### Description

---

This Anti-FOXP3 [ZN248] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                               |
|-----------------------|-----------------------------------------------|
| <b>SKU:</b>           | AGMB06293                                     |
| <b>Contents:</b>      | 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml) |
| <b>Synonyms:</b>      | JM2, AIID, IPEX, PIDX, XPID, DIETER           |
| <b>Applications:</b>  | <b>IHC-P</b>                                  |
| <b>Research Area:</b> | Cell Biology                                  |
| <b>Form:</b>          | Liquid                                        |

### Antibody Data

---

|                      |                                                                                         |
|----------------------|-----------------------------------------------------------------------------------------|
| <b>Reactivity:</b>   | Human                                                                                   |
| <b>Clone:</b>        | ZN248                                                                                   |
| <b>Clonality:</b>    | Monoclonal Antibody                                                                     |
| <b>SwissProt ID:</b> | Q9BZS1                                                                                  |
| <b>Immunogen:</b>    | Full length human recombinant protein of human FOXP3 (NP_001107849) produced in E.coli. |
| <b>Gene ID:</b>      | 50943                                                                                   |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** FOXP3

**Host Species:** Mouse

**Isotype:** IgG2a

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | IHC                | 1:50-1:100                     |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.