

Anti-Cytokeratin 7 [R03-4W-4] Monoclonal Antibody

AGMB06366

Description

This Anti-Cytokeratin 7 [R03-4W-4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06366
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CK 7, CK-7, ck7, Cytokeratin 7, Cytokeratin-7, Cytokeratin7, D15Wsu77e, K2C7, K2C7_HUMAN, K7, Keratin 55k type ii cytoskeletal, Keratin 7, Keratin simple epithelial type 1 k7, Keratin type II cytoskeletal 7
Applications:	IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-4W-4
Clonality:	Monoclonal Antibody
SwissProt ID:	P08729
Immunogen:	A synthesized peptide derived from human KRT7
Gene ID:	3855

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: KRT7
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.