

Anti-IL-2 Receptor alpha [R03-5O-2] Monoclonal Antibody

AGMB06402

Description

This Anti-IL-2 Receptor alpha [R03-5O-2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06402
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	IL2RA, Interleukin-2 receptor subunit alpha, IL-2 receptor subunit alpha, IL-2-RA, IL-2R subunit alpha, IL2-RA, TAC antigen, p55, CD25
Applications:	IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-5O-2
Clonality:	Monoclonal Antibody
SwissProt ID:	P01589
Immunogen:	A synthesized peptide derived from human IL2RA
Gene ID:	3559

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: IL2RA
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.