

Anti-Human IgA [R01-8W-2] Monoclonal Antibody

AGMB06439

Description

This Anti-Human IgA [R01-8W-2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

| | |
|-----------------------|---|
| SKU: | AGMB06439 |
| Contents: | 50µl Bradford Reagent: 1 vial (2ml) |
| Synonyms: | IGHA1, Ig alpha-1 chain C region, IGH A2, Ig alpha-2 chain C region |
| Applications: | IHC-P |
| Research Area: | Immunology |
| Form: | Liquid |

Antibody Data

| | |
|----------------------|--|
| Reactivity: | Human |
| Clone: | R01-8W-2 |
| Clonality: | Monoclonal Antibody |
| SwissProt ID: | P01876/P01877 |
| Immunogen: | A synthesized peptide derived from human IGH A1/IGH A2 |
| Gene ID: | 3493/3494 |
| Gene Name: | IGH A1/IGH A2 |

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|--------------------------|-------------|-------------------------|
| | IHC | 1:100-1:200 |

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.