

# Anti-BRG1 [R01-3P-8] Monoclonal Antibody

## AGMB06442

### Description

---

This Anti-BRG1 [R01-3P-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB06442
<b>Contents:</b>	50µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	SMARCA4, BAF190A, BRG1, SNF2B, SNF2L4, Transcription activator BRG1, ATP-dependent helicase SMARCA4, BRG1-associated factor 190A, BAF190A, Mitotic growth and transcription activator, Protein BRG-1, Protein brahma homolog 1, SNF2-beta, SWI/S
<b>Applications:</b>	<b>IHC-P</b>
<b>Research Area:</b>	Neuroscience
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R01-3P-8
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P51532
<b>Immunogen:</b>	A synthesized peptide derived from human SMARCA4

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 6597  
**Gene Name:** SMARCA4  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.