

Anti-TIM3 [R01-9S-5] Monoclonal Antibody

AGMB06451

Description

This Anti-TIM3 [R01-9S-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06451
Contents:	50µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD366, KIM-3, SPTCL, TIMD3, Tim-3, TIMD-3, HAVcr-2, TIM-3/HAVCR2
Applications:	IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R01-9S-5
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8TDQ0
Immunogen:	A synthesized peptide derived from human HAVCR2
Gene ID:	84868
Gene Name:	HAVCR2

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.