

Anti-LAG 3 [R05-5R-3] Monoclonal Antibody

AGMB06452

Description

This Anti-LAG 3 [R05-5R-3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06452
Contents:	50µl Bradford Reagent: 1 vial (2ml)
Synonyms:	LAG3, FDC, Lymphocyte activation gene 3 protein, LAG-3, Protein FDC, CD223
Applications:	IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-5R-3
Clonality:	Monoclonal Antibody
SwissProt ID:	P18627
Immunogen:	A synthesized peptide derived from human LAG3
Gene ID:	3902

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: LAG3
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.