

# Anti-AR-V7 [R09-2E2] Monoclonal Antibody

AGMB06463

## Description

---

This Anti-AR-V7 [R09-2E2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

|                       |  |
|-----------------------|--|
| <b>SKU:</b>           | AGMB06463  |
| <b>Contents:</b>      | 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)  |
| <b>Synonyms:</b>      | AR,AR3,DHTR, NR3C4, Dihydrotestosterone receptor, Nuclear receptor subfamily 3 group C member 4,Androgen receptor splice variant 3 |
| <b>Applications:</b>  | <b>IHC-P</b>   |
| <b>Research Area:</b> | Epigenetics  |
| <b>Form:</b>          | Liquid   |

## Antibody Data

---

|                      |   |
|----------------------|---|
| <b>Reactivity:</b>   | Human   |
| <b>Clone:</b>        | R09-2E2   |
| <b>Clonality:</b>    | Monoclonal Antibody                             |
| <b>SwissProt ID:</b> | C0JKD3  |
| <b>Immunogen:</b>    | A synthesized peptide derived from human AR-V7. |
| <b>Gene ID:</b>      | 367   |

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** AR  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | IHC                | 1:100-1:200                    |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.