

Anti-CD23 [R04-3O-1] Monoclonal Antibody

AGMB06499

Description

This Anti-CD23 [R04-3O-1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06499
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FCER2, CD23A, CLEC4J, FCE2, IGEBF, Low affinity immunoglobulin epsilon Fc receptor, BLAST-2, C-type lectin domain family 4 member J, Fc-epsilon-R1I, Immunoglobulin E-binding factor, Lymphocyte IgE receptor, CD23
Applications:	IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R04-3O-1
Clonality:	Monoclonal Antibody
SwissProt ID:	P06734
Immunogen:	Recombinant protein of human CD23
Gene ID:	2208

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: FCER2
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	IHC-P	1:500-1:2000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.