

Anti-alpha Smooth Muscle Actin [R09-7F-5] Monoclonal Antibody

AGMB06520

Description

This Anti-alpha Smooth Muscle Actin [R09-7F-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06520
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ACTA2, ACTSA, ACTVS, GIG46, Actin, aortic smooth muscle, Alpha-actin-2, Cell growth-inhibiting gene 46 protein, α-SMA
Applications:	IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-7F-5
Clonality:	Monoclonal Antibody
SwissProt ID:	P62736
Immunogen:	A synthesized peptide derived from human alpha smooth muscle Actin
Gene ID:	59

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ACTA2
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:500-1:2000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.