

Anti-CD3 [R05-7H-4] Monoclonal Antibody

AGMB06527

Description

This Anti-CD3 [R05-7H-4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06527
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD3D, CD3E, CD3G, CD3Z,
Applications:	IHC-P
Research Area:	-
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-7H-4
Clonality:	Monoclonal Antibody
SwissProt ID:	P04234/P07766/P09693/P20963
Immunogen:	A synthesized peptide derived from human CD3
Gene ID:	915, 916, 917, 919
Gene Name:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:200-1:2000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.