

# Anti-FOXP3 [8R9-I5-H5] Monoclonal Antibody

## AGMB06531

### Description

---

This Anti-FOXP3 [8R9-I5-H5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB06531
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	JM2, AIID, IPEX, PIDX, XPID, DIETER
<b>Applications:</b>	<b>IHC-P</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	8R9-I5-H5
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q9BZS1
<b>Immunogen:</b>	-
<b>Gene ID:</b>	50943
<b>Gene Name:</b>	FOXP3

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Purification  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:200-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.