

Anti-PVRL4 [1H8-K5-I9] Monoclonal Antibody

AGMB06532

Description

This Anti-PVRL4 [1H8-K5-I9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06532
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	LNIR, PRR4, EDSS1, nectin-4
Applications:	WB IHC-P
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	1H8-K5-I9
Clonality:	Monoclonal Antibody
SwissProt ID:	Q96NY8
Immunogen:	-
Gene ID:	81607
Gene Name:	NECTIN4

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purification

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 56 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:200-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.