

Anti-VEGFA [6E9-W3-A4] Monoclonal Antibody

AGMB06538

Description

This Anti-VEGFA [6E9-W3-A4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06538
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	VEGFA, VEGF, Vascular endothelial growth factor A, VEGF-A, Vascular permeability factor, VPF, Folliculostellate cell-derived growth factor, Glioma-derived endothelial cell mitogen, MGC70609,
Applications:	IHC-P
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	6E9-W3-A4
Clonality:	Monoclonal Antibody
SwissProt ID:	P15692
Immunogen:	-
Gene ID:	7422

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: VEGFA
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purification
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.