

Anti-FBXO11 Monoclonal Antibody -Knockout Validated

AGMB06563

Description

This Anti-FBXO11 Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06563
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FBX11, PRMT9, UBR6, UG063H01, VIT1
Applications:	WB
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	1000
Clonality:	Monoclonal Antibody
SwissProt ID:	Q86XK2
Immunogen:	Purified recombinant human FBXO11 protein fragments expressed in E.coli.
Gene ID:	80204

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: FBXO11

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 104 kDa, Observed MW: 104 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 1% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.