

Anti-E2F4 [R09-7E6] Monoclonal Antibody -Knockout Validated

AGMB06566

Description

This Anti-E2F4 [R09-7E6] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06566
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	E2F4, Transcription factor E2F4, E2F-4, E2F 4, E2F transcription factor 4, E2F transcription factor 4 p107/p130 binding,
Applications:	WB IHC-P
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-7E6
Clonality:	Monoclonal Antibody
SwissProt ID:	Q16254
Immunogen:	A synthetic peptide of human E2F4
Gene ID:	1874

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: E2F4
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 44 kDa, Observed MW: 65 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:1000-1:5000
	IHC-P	1:200-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.