

Anti-p150 CAF1 [R06-3A5] Monoclonal Antibody -Knockout Validated

AGMB06582

Description

This Anti-p150 CAF1 [R06-3A5] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

| | |
|-----------------------|---|
| SKU: | AGMB06582 |
| Contents: | 20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml) |
| Synonyms: | CHAF1A, CAF, CAF1P150, Chromatin assembly factor 1 subunit A, CAF-1 subunit A, Chromatin assembly factor I p150 subunit, CAF-I 150 kDa subunit, CAF-I p150, hp150 |
| Applications: | WB IHC-P ICC/IF FC |
| Research Area: | Epigenetics |
| Form: | Liquid |

Antibody Data

| | |
|----------------------|--|
| Reactivity: | Human |
| Clone: | R06-3A5 |
| Clonality: | Monoclonal Antibody |
| SwissProt ID: | Q13111 |
| Immunogen: | A synthetic peptide of human p150 CAF1 |
| Gene ID: | 10036 |

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CHAF1A
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 107 kDa, Observed MW: 150 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|---------------------------------|--------------------|--------------------------------|
| | WB | 1:1000-1:5000 |
| | IHC-P | 1:100-1:200 |
| | ICC/IF | 1:20-1:50 |
| | FC | 1:20-1:50 |

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.