

Anti-YAP1 [6B11] Monoclonal Antibody -Knockout Validated

AGMB06607

Description

This Anti-YAP1 [6B11] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06607
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	YAP1, YAP65, Yorkie homolog, 65 kDa Yes-associated protein, YAP65
Applications:	WB
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	6B11
Clonality:	Monoclonal Antibody
SwissProt ID:	P46937
Immunogen:	Purified recombinant human YAP1 protein fragments expressed in E.coli.
Gene ID:	10413
Gene Name:	YAP1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 54 kDa, Observed MW: 70 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 1% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.