

Anti-GIPC1 [5C10] Monoclonal Antibody -Knockout Validated

AGMB06651

Description

This Anti-GIPC1 [5C10] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06651
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	NIP, GIPC, IIP-1, OPDM2, TIP-2, SEMCAP, C19orf3, Hs.6454, GLUT1CBP, RGS19IP1, SYNECTIN, SYNECTIIN, GIPC1
Applications:	WB IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	5C10
Clonality:	Monoclonal Antibody
SwissProt ID:	O14908
Immunogen:	Full length human recombinant protein of human GIPC1 produced in E.coli.
Gene ID:	10755

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: GIPC1
Host Species: Mouse
Isotype: IgG2a
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 36 kDa, Observed MW: 42 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 1% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.