

Anti-L1CAM [1H3] Monoclonal Antibody -Knockout Validated

AGMB06699

Description

This Anti-L1CAM [1H3] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06699
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	L1CAM, CAML1, MIC5, Neural cell adhesion molecule L1, N-CAM-L1, NCAM-L1, CD antigen CD171
Applications:	WB FC
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	1H3
Clonality:	Monoclonal Antibody
SwissProt ID:	P32004
Immunogen:	Full length human recombinant protein of human L1CAM produced in 293T cell.
Gene ID:	3897

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: L1CAM
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 140 kDa, Observed MW: 250 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 1% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.