

Anti-PP-T [6F1] Monoclonal Antibody -Knockout Validated

AGMB06723

Description

This Anti-PP-T [6F1] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06723
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	PPP5C, FLJ36922, FLJ55954, OTTHUMP00000165834, PP-T, PP5, PPP5, PPP5_HUMAN, PPT, Protein phosphatase 5, catalytic subunit, Protein phosphatase T, Serine/threonine protein phosphatase 5, Serine/threonine-protein phosphatase 5.
Applications:	WB IHC-P
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	6F1
Clonality:	Monoclonal Antibody
SwissProt ID:	P53041
Immunogen:	Full length human recombinant protein of human PPP5C produced in 293T cell.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 5536

Gene Name: PPP5C

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 57 kDa, Observed MW: 57 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 1% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.