

# Anti-STAT3 [210000000] Monoclonal Antibody -Knockout Validated

## AGMB06745

### Description

---

This Anti-STAT3 [210000000] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB06745
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	STAT3, APRF, Signal transducer and activator of transcription 3, Acute-phase response factor, 1110034C02Rik, Acute Phase Response Factor, Acute-phase response factor,
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b>
<b>Research Area:</b>	Signal Transduction
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	210000000
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P40763
<b>Immunogen:</b>	Purified recombinant human STAT3 protein fragments expressed in E.coli.
<b>Gene ID:</b>	6774

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** STAT3

**Host Species:** Mouse

**Isotype:** IgG1

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 88 kDa, Observed MW: 88 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 1% BSA and 0.02% sodium azide, pH 7.3.

	Application	Antibody Dilution Ratio
<b>Antibody Dilution Ratio:</b>	WB	1:500-1:1000
	IHC	1:50-1:100

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.