

Anti-Integrin alpha V [R02-6G3] Monoclonal Antibody -Knockout Validated

AGMB06771

Description

This Anti-Integrin alpha V [R02-6G3] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06771
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ITGAV, MSK8, VNRA, Integrin alpha-V, Vitronectin receptor subunit alpha, CD antigen CD51
Applications:	WB IHC-P ICC/IF IP
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R02-6G3
Clonality:	Monoclonal Antibody
SwissProt ID:	P06756
Immunogen:	Recombinant protein of human Integrin alpha V
Gene ID:	3685

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ITGAV
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 116 kDa, Observed MW: 130 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:200-1:500
	ICC/IF	1:200-1:500
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.