

# Anti-Phospho-Presenilin 1 (Ser310) [R07-6E3] Monoclonal Antibody - Knockout Validated

AGMB06773

## Description

---

This Anti-Phospho-Presenilin 1 (Ser310) [R07-6E3] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB06773
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	AD3, FAD, PS1, PS-1, S182, PSNL1, ACNINV3
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Neuroscience
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R07-6E3
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P49768
<b>Immunogen:</b>	A synthetic phosphopeptide corresponding to residues surrounding Ser310 of human Presenilin 1
<b>Gene ID:</b>	5663

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** PSEN1  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Phosphorylated  
**Molecular Weight:** Calculated MW: 53 kDa, Observed MW: 22 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

	Application	Antibody Dilution Ratio
<b>Antibody Dilution Ratio:</b>	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.