

Flag Tag Affinity Gel

AGMB06809

Description

This Flag Tag Affinity Gel is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06809
Contents:	2mL, 5mL Bradford Reagent: 1 vial (2ml)
Synonyms:	DYKDDDDK Tag, Flag, FLAG
Applications:	IP Purification
Research Area:	Tags & Cell Markers
Form:	Slurry

Antibody Data

Reactivity:	Species-Independent
Clone:	-
Clonality:	Monoclonal Antibody
SwissProt ID:	-
Immunogen:	DYKDDDDK
Gene ID:	-
Gene Name:	-
Host Species:	Mouse

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Isotype: IgG2b
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at -20°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: 50% (v/v) in PBS, pH 7.4, containing 0.02% Sodium Azide and 50% glycerol.

Antibody Dilution Ratio: -

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.