

# Anti-PRMT1 (4E9) [4E9-H9-H12] Monoclonal Antibody

## AGMB06819

### Description

---

This Anti-PRMT1 (4E9) [4E9-H9-H12] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB06819
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	ANM 1, ANM1, ANM1_HUMAN, complete cds, HCP 1, HCP1, Heterogeneous nuclear ribonucleoprotein methyltransferase 1 like 2, Heterogeneous nuclear ribonucleoproteins methyltransferase like 2, Heterogeneous nuclear ribonucleoproteins methyltransferase like2, Histone-arginine N-methyltransferase PRMT1, HMT 2, HMT1 (hnRNP methyltransferase, HMT1 (hnRNP methyltransferase S. cerevisiae) like 2, HMT1 hnRNP methyltransferase
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	4E9-H9-H12
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q99873

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Immunogen:** Purified recombinant human PRMT1 protein fragments expressed in E.coli.

**Gene ID:** 3276

**Gene Name:** PRMT1

**Host Species:** Mouse

**Isotype:** IgG2b

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 42 kDa, Observed MW: 42 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.