

# Anti-beta Actin [8F10-G10] Monoclonal Antibody (HRP)

AGMB06827

## Description

---

This Anti-beta Actin [8F10-G10] Monoclonal Antibody (HRP) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB06827
<b>Contents:</b>	1mL, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	ACTIN, ACTB, BETA-ACTIN, beta actin, actin beta
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Isotype/Loading Controls
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat, Hamster, Monkey, Goat, Rabbit, Zebrafish, Saccharomycescerevisiae, Arabidopsis, Pichiapastoris, Drosophila, LeafofRice, FruitofCucumber, ChlamydomonasReinhardtii
<b>Clone:</b>	8F10-G10
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P60709
<b>Immunogen:</b>	Purified recombinant human beta-actin protein fragments expressed in E.coli.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 60

**Gene Name:** ACTB

**Host Species:** Mouse

**Isotype:** IgG2b

**Purification:** Affinity Purified

**Conjugated:** HRP

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 42 kDa, Observed MW: 42 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:5000-1:10000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.