

Anti-GFP Tag [R04-6A1] Monoclonal Antibody

AGMB06832

Description

This Anti-GFP Tag [R04-6A1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06832
Contents:	1mL, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	GFP, Green Fluorescent Protein, enhanced Green Fluorescent Protein.
Applications:	WB IHC-P ICC/IF FC
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Species-Independent
Clone:	R04-6A1
Clonality:	Monoclonal Antibody
SwissProt ID:	P42212
Immunogen:	Recombinant protein of GFP
Gene ID:	-
Gene Name:	GFP

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 27 kDa, Observed MW: Refer to figures

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC	1:100-1:200
	ICC/IF	1:200-1:500
	FC	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.