

Anti-c-Myc Tag [9F6-7B4-4D8] Monoclonal Antibody

AGMB06848

Description

This Anti-c-Myc Tag [9F6-7B4-4D8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06848
Contents:	1mL, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	c-Myc, Myc tag, Myc protein, Myc epitope tag, Myc2
Applications:	WB
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Species-Independent
Clone:	9F6-7B4-4D8
Clonality:	Monoclonal Antibody
SwissProt ID:	-
Immunogen:	Synthetic peptide EQKLISEEDL coupled to KLH.
Gene ID:	-
Gene Name:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
		WB

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.