

# Anti-c-Myc Tag [7D9-4C6-6F10] Monoclonal Antibody

## AGMB06855

### Description

---

This Anti-c-Myc Tag [7D9-4C6-6F10] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB06855
<b>Contents:</b>	1mL, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	c-Myc, Myc tag, Myc protein, Myc epitope tag, Myc2
<b>Applications:</b>	<b>WB</b> <b>ICC/IF</b> <b>IP</b> <b>ELISA</b>
<b>Research Area:</b>	Tags & Cell Markers
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Species-Independent
<b>Clone:</b>	7D9-4C6-6F10
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	-
<b>Immunogen:</b>	Synthetic Peptide of Myc-Tag
<b>Gene ID:</b>	-
<b>Gene Name:</b>	-

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse  
**Isotype:** IgG1  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:1000-1:5000
	IP	1-1:200
	IF	1:1000
	ELISA	1:5000-1:20000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.