

Anti-GFP Tag (9A7) [9A7-1E6-4G1] Monoclonal Antibody

AGMB06858

Description

This Anti-GFP Tag (9A7) [9A7-1E6-4G1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06858
Contents:	1mL, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	GFP, Green Fluorescent Protein, enhanced Green Fluorescent Protein.
Applications:	WB
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Species-Independent
Clone:	9A7-1E6-4G1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q963I9
Immunogen:	Recombinant Protein of GFP-Tag
Gene ID:	-
Gene Name:	GFP

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 27 kDa, Observed MW: Refer to figures

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.