

Anti-Flag Tag (10E1) [10E1-1D8-5F5] Monoclonal Antibody

AGMB06871

Description

This Anti-Flag Tag (10E1) [10E1-1D8-5F5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06871
Contents:	1mL, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	DYKDDDDK Tag, Flag, FLAG
Applications:	WB IP
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Species-Independent
Clone:	10E1-1D8-5F5
Clonality:	Monoclonal Antibody
SwissProt ID:	-
Immunogen:	Monoclonal antibodies are produced by immunizing mouse with synthetic peptide (DYKDDDDK), coupled to KLH.
Gene ID:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: -
Host Species: Mouse
Isotype: IgG2b
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:10000
	IP	1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.