

# Anti-HA Tag Monoclonal Antibody (HRP)

AGMB06876

## Description

---

This Anti-HA Tag Monoclonal Antibody (HRP) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB06876
<b>Contents:</b>	1mL, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	HA epitope tag, HA1, HA2, hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain
<b>Applications:</b>	<b>WB</b> <b>ELISA</b>
<b>Research Area:</b>	Tags & Cell Markers
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Species-Independent
<b>Clone:</b>	-
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	-
<b>Immunogen:</b>	YPYDVPDYA conjugated to KLH.
<b>Gene ID:</b>	-

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** -  
**Host Species:** Mouse  
**Isotype:** IgG1  
**Purification:** Affinity Purified  
**Conjugated:** HRP  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:1000-1:10000
	ELISA	1:5000-1:20000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.