

Anti-PCSK9 Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB06901

Description

This Anti-PCSK9 Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06901
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	NARC-1,PC9,Subtilisin/kexin-like protease PC9,Narc1
Applications:	ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	-
Clonality:	Monoclonal Antibody
SwissProt ID:	Q80W65
Immunogen:	-
Gene ID:	100102
Gene Name:	Pcsk9

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit
Isotype: -
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 2-8°C. Avoid freezing.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.