

Anti-HA tag [R04-3B-1] Monoclonal Antibody

AGMB06905

Description

This Anti-HA tag [R04-3B-1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06905
Contents:	1mL, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	HA epitope tag, HA1, HA2, hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Species-Independent
Clone:	R04-3B-1
Clonality:	Monoclonal Antibody
SwissProt ID:	-
Immunogen:	HA synthetic peptide conjugated to KLH.
Gene ID:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: -
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:1000-1:5000
	ICC/IF	1:500-1:1000
	FC	1:200-1:500
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.